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Table 1. Clinical breakpoints and screening cut-off values for carbapenemase-producing Enterobacteriaceae (according to EUCAST methodology).

Carbapenem	MIC (mg/L)		Disk diffusion zone diameter (mm) with 10 µg disks	
	S/I breakpoint	Screening cut-off	S/I breakpoint	Screening cut-off
Meropenem ¹	≤2	>0.125	≥22	<28 ²
Ertapenem ³	≤0.5	>0.125	≥25	<25

¹Best balance of sensitivity and specificity

²Isolates with 25-27 mm only need to be investigated for carbapenemase-production if they are resistant to piperacillin-tazobactam and/or temocillin (temocillin contributes more to the specificity). Investigation for carbapenemases is always warranted if zone diameter of meropenem is <25 mm.

³High sensitivity but low specificity. Can be used as an alternative screening agent, but isolates with ESBL and AmpC may be resistant without having carbapenemases.

Kinnitav test

T1. AmpC kinnitav test (tsefsalosporiinid +/-kloksatsilliin).

T2. ESBL_A test. Tseftasidiimi ja tsefotaksiimi klavulaanhappe sünergismi test.

Method	Antimicrobial agent (disk content)	ESBL confirmation is positive if
ESBL gradient test	Cefotaxime +/- clavulanic acid	MIC ratio ≥ 8 or deformed ellipse present
	Ceftazidime +/- clavulanic acid	MIC ratio ≥ 8 or deformed ellipse present
Combination disk diffusion test (CDT)	Cefotaxime (30 μ g) +/- clavulanic acid (10 μ g)	≥ 5 mm increase in inhibition zone
	Ceftazidime (30 μ g) +/- clavulanic acid (10 μ g)	≥ 5 mm increase in inhibition zone
Broth microdilution	Cefotaxime +/- clavulanic acid (4 mg/L)	MIC ratio ≥ 8
	Ceftazidime +/- clavulanic acid (4 mg/L)	MIC ratio ≥ 8
	Cefepime +/- clavulanic acid (4 mg/L)	MIC ratio ≥ 8
Double disk synergy test (DDST)	Cefotaxime, ceftazidime and cefepime	Expansion of indicator cephalosporin inhibition zone towards amoxicillin-clavulanic acid disk

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Method	Antibiotic	Confirmation is positive if
ESBL gradient test Etest [®] ESBL	Cefepime +/- clavulanic acid	MIC ratio ≥ 8 or deformed ellipse present
Combination disk diffusion test	Cefepime (30 μ g) +/- clavulanic acid (10 μ g)	≥ 5 mm increase in inhibition zone
Broth microdilution	Cefepime +/- clavulanic acid (fixed concentration 4 mg/L)	MIC ratio ≥ 8
Double disk synergy test (DDST)	Cefotaxime, ceftazidime, Cefepime	Expansion of indicator cephalosporin inhibition zone towards amoxicillin-clavulanic acid disk

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¹ Temocillin susceptibility test is recommended only in cases where no synergy is detected, in order to differentiate between ESBL + porin loss and OXA-48-like enzymes (23, 24). When other enzymes are present the susceptibility is variable and does not provide any further indication of the β-lactamase present.

² There is one report supporting the use of commercial tablets containing double inhibitors (DPA or EDTA plus APBA or PBA) (25), but multi-centre studies or multiple single centre studies are lacking. This combination confers high-level resistance to carbapenems and is rare outside Greece.

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