



# EUCAST Rapid Antimicrobial Susceptibility Testing (RAST) directly from positive blood cultures

Emma Jonasson, Biomedical scientist  
EUCAST development laboratory  
Växjö, Sweden

Tallinn, Estonia 4-6 September 2024

# RAST methodology and performance

# Why a standardised rapid method?

- A rapid (and correct) AST result is crucial in the treatment of septicaemia.
- Many microbiology laboratories have developed in-house methods for rapid AST directly from positive blood cultures.
- EUCAST wished to take responsibility for validating **one** RAST method for positive blood cultures to harmonize between laboratories and ensure reliable results.

# EUCAST RAST validation

- Spiked blood culture bottles
  - Isolates with various levels of susceptibility/resistance
  - Tests performed at EDL
  - MIC with BMD was used as a reference
- Multi-laboratory trial
  - Local clinical isolates
  - Local BC systems, MH media and disks
  - Northern Europe 2017: 40 laboratories
  - Southern Europe 2018: 15 laboratories

# EUCAST RAST methodology

- EUCAST standard disk diffusion method with modifications of:
  - Inoculum: 100-150  $\mu$ L, positive blood culture
  - Incubation time: 4, 6 and/or 8 h (depending on species)
- Developed for:
  - *E. coli*, *K. pneumoniae*, *S. enterica*, *P. aeruginosa*, *A. baumannii*,  
*S. aureus*, *E. faecalis*, *E. faecium* and *S. pneumoniae*

# EUCAST RAST methodology

- EUCAST standard disk diffusion method with modifications of:
  - Inoculum: 100-150  $\mu$ L, positive blood culture
  - Incubation time: 4, 6 and/or 8 h (depending on species) **and 16-20 h**
- Developed for:
  - *E. coli*, *K. pneumoniae*, *S. enterica*, *P. aeruginosa*, *A. baumannii*,  
*S. aureus*, *E. faecalis*, *E. faecium* and *S. pneumoniae*

# Incubation

- 4, 6 and/or 8 hours
  - Read plates within  $\pm 5$  minutes of the specified incubation time
  - Re-incubate the plate within 10 minutes
- 16-20 hours intervall
  - Use when it is not possible to read results after 4, 6 and/or 8 hours
  - Breakpoints available for all species


# Reading of plates

- General reading instructions
  - Hold the plate at 30 cm distance
- 4, 6 and 8 hours
  - Read plates from the front without lid
  - Read zones only when a clear zone edge is visible
- 16-20 hours
  - Read MH plates from the back of the plate
  - Read MH-F plates from the front without lid




# Reading of zones

4, 6 and/or 8 hours



*E. coli* at 4h incubation

16-20 hours



*E. coli* at 20h incubation

# Readability (%)

The proportion of zone diameters (%) which are possible to read after 4 – 20 h of incubation.

Organism	4 hours (%)	6 hours (%)	8 hours (%)	16-20 hours
<i>Escherichia coli</i>	90	99	99	100
<i>Klebsiella pneumoniae</i>	96	98	98	100
<i>Salmonella enterica</i>	93	100	100	100
<i>Pseudomonas aeruginosa</i>	-	88	97	100
<i>Acinetobacter baumannii</i>	99	100	100	ND
<i>Staphylococcus aureus</i>	55	91	95	100
<i>Enterococcus faecalis</i>	93	99	100	ND
<i>Enterococcus faecium</i>	44	93	99	ND
<i>Streptococcus pneumoniae</i>	68	83	95	100

# Readability (%)

The proportion of zone diameters (%) which are possible to read after 4 – 20 h of incubation.

Organism	4 hours (%)	6 hours (%)	8 hours (%)	16-20 hours
<i>Escherichia coli</i>	90	99	99	100
<i>Klebsiella pneumoniae</i>	96	98	98	100
<i>Salmonella enterica</i>	93	100	100	100
<i>Pseudomonas aeruginosa</i>	-	88	97	100
<i>Acinetobacter baumannii</i>	99	100	100	ND
<i>Staphylococcus aureus</i>	55	91	95	100
<i>Enterococcus faecalis</i>	93	99	100	ND
<i>Enterococcus faecium</i>	44	93	99	ND
<i>Streptococcus pneumoniae</i>	68	83	95	100

# Readability (%)

The proportion of zone diameters (%) which are possible to read after 4 – 20 h of incubation.







Organism	4 hours (%)	6 hours (%)	8 hours (%)	16-20 hours
<i>Escherichia coli</i>	90	99	99	100
<i>Klebsiella pneumoniae</i>	96	98	98	100
<i>Salmonella enterica</i>	93	100	100	100
<i>Pseudomonas aeruginosa</i>	-	88	97	100
<i>Acinetobacter baumannii</i>	99	100	100	ND
<i>Staphylococcus aureus</i>	55	91	95	100
<i>Enterococcus faecalis</i>	93	99	100	ND
<i>Enterococcus faecium</i>	44	93	99	ND
<i>Streptococcus pneumoniae</i>	68	83	95	100

# Interpretation of results







- Interpret inhibition zone diameters according to the latest version of the RAST breakpoint tables.
- Leave the report blank if:
  - Cannot read the zone in a reliable way
  - The zone diameter is in the ATU.

Only interpret results for species with RAST breakpoints!







# *E. coli* with cefotaxime 5 µg vs. BMD



# *E. coli* with cefotaxime 5 µg vs. BMD




# *E. coli* with cefotaxime 5 µg vs. BMD





# Area of Technical Uncertainty (ATU)



- No AST report was given for results within ATU
  - ATU after 4 or 6 hours = plates are reincubated.
  - ATU also after 8 hours, retest with standard methodology.
  - ATU after 16-20 hours, retest with standard methodology.

# Categorical agreement

Northern and Southern European study, clinical samples and all species.

Number of readable zones 4h: 5811, 6h: 6921, 8h: 6561

Breakpoint version 1.

Incubation time		4h	6h	8h
Results calculated on readable zones (%)				
ATU		16	7.5	5.7
Interpreted o S or R		84	93	94
Errors calculated on zones interpreted to S or R (%)				
	mE	0.6	0.6	0.8
Errors	ME	2.1	1.1	0.9
	VME	0.2	0.4	0.5
		3.0	2.1	2.2

# Quality Controls (QC)

1. Standard disk diffusion QC – to control AST materials
2. RAST QC - for implementation and changes of the RAST method
  1. From BC bottles with added blood
  2. RAST methodology 4, 6, 8 and 16-20 hours reading
  3. Practice reading

RAST in a clinical laboratory

# Implementation

- Should RAST be performed on all bottles?
- Which incubation time?
- How is the results interpreted?
  - Avoid manual interpretation
  - Incorporate in LIS system
- Staff training
  - Theoretical
  - Practical
- Information to clinician



Incubation start	4h reading	6h reading
8.00	12.00	14.00
9.00	13.00	15.00

# Daily routine

- Priorities - allocate staff resources
  - Much work in parallel (Gram stain, MALDI-TOF, plating, phone, RAST)
- Read at correct time
  - Timer
  - Whiteboard




LABNUMBER	TID	4h/6h
563335 g-	8 <sup>55</sup>	12 <sup>55</sup>
563780 staf	9 <sup>30</sup>	13 <sup>30</sup>

# Results from our clinical lab

- Primary reading time 4h
  - First AST results in 5.5 hours from removal of bottle
- Fraction of isolates with completed RAST, in average 78 %

# Fraction of isolates with completed RAST

Time period 14 months, RAST 4, 6 and 8 hours  
N=1300 isolates









# Categorical agreement


RAST 4, 6 and 8 hours during 1 year

Species		<i>E. coli</i>			<i>S. aureus</i>		
		4h	6h	8h	4h	6h	8h
Incubation time		4h	6h	8h	4h	6h	8h
Number of isolates tested		138	223	4	68	139	3
Readable zones (%)		98	100	75	95	98	100
Results calculated on readable zones (%)							
Not interpreted to S or R (ATU)		17	11	10	8	10	0
Interpreted to S		79	83	90	91	87	100
Interpreted to R		4	6	0	1	4	0
Errors calculated on zones interpreted to S or R (%)							
Errors	mE	1.3	0.9	0	0	0	0
	ME	0.1	0.5	0	0	0.8	0
	VME	0	0	0	0	0	0
	Total errors:	1.4	1.4	0	0	0.8	0

# EUCAST website

← ↻ 🏠 🔒 https://www.eucast.org/rapid\_ast\_in\_bloodcultures 🔍 🔍 🔍 ⚙️ 📄 ⭐ 📧 📱 📺

Home Contact Sitemap Newsletter    

 **EUCAST** EUROPEAN COMMITTEE ON ANTIMICROBIAL SUSCEPTIBILITY TESTING  
European Society of Clinical Microbiology and Infectious Diseases

search term 🔍

## Rapid AST in bloodcultures

Organization

Public consultations

EUCAST News

Definitions of S, I and R

Clinical breakpoints and dosing

**Rapid AST in blood cultures**

Calibration files

RAST from EUCAST in publications


Methods

QC

Breakpoints for short incubation

Screening for resistance mechanisms

FAQ on RAST



Rapid AST in bloodcultures ▾

### Rapid AST directly from blood culture bottles

April 12, 2022

EUCAST has developed a method for rapid AST (reading at 4, 6 or 8h and since April 2022 also after 16-20 hours incubation) directly from positive blood culture bottles (RAST). The [rationale is available in JAC](#).

Following the initial development, [published in 2019](#), a clinical trial in 55 laboratories was performed and [published in 2020](#).

These are the essential steps in the RAST method:

# Future

- Breakpoints for new agents and/or breakpoints for additional species when considered relevant by the EUCAST committee
- A similar methodology for urinary tract samples from patients with bloodstream infections – to reduce the time to susceptibility results

# Questions?



[Emma.k.jonasson@kronoberg.se](mailto:Emma.k.jonasson@kronoberg.se)