

# Point of Care Cardiac Troponin I and T Assay Analytical Characteristics Designated by Manufacturer

## IFCC Committee on Clinical Applications of Cardiac Bio-Markers (C-CB) v062024

Company/Platform/ Assay	LoB (µg/L unless noted)	LoD (µg/L unless noted)	% CV at 99 <sup>th</sup> Percentile	Conc at 20% CV (µg/L unless noted)	Conc at 10% CV (µg/L unless noted)	Reference Population N, Ages, Sex	Specimen Type	Calibrators used	99 <sup>th</sup> Percentile (µg/L unless noted)	Percent Normals Measured ≥ LoD	Statistic Used to Calc 99 <sup>th</sup> Percentile	Epitopes Recognized by Antibodies	Country of Package Insert: Version Date
Abbott i-STAT	0.02	NP	16.5%	0.07	0.1	Overall n = 162	Sodium and lithium heparinized whole blood and plasma	NP	Overall: 0.08	NP	NP	NP	US: Rev. Date: 01-Jul-13
LSI Medience PATHFAST cTnI; commercial	NP	1 ng/L	< 6%	2 ng/L	3.1 ng/L	Overall n =474 18-86y  F: 236 M: 238	Heparin-Na, heparin-Li or EDTA whole blood or plasma	Human cardiac I-T- C complex	Overall: 15.48 ng/L  F: 11.46 ng/L M: 16.91 ng/L	Overall: 76.3%	Non-Parametric	C: 41-49, D: 71-116, 163-209	WW except US & Japan: Ver.6, Apr. 2021
LSI Medience PATHFAST hs-cTnI /PATHFAST cTnI-II	1.23 ng/L	2.33 ng/L	6.1%	4 ng/L	15 ng/L	Overall n=734 Age >18  F: 352 M: 382	Heparin-Na, heparin-Li or EDTA whole blood or plasma	Human cardiac I-T- C complex	Overall: 27.9 ng/L  F: 20.3 ng/L M: 29.7 ng/L	Overall: 66.3%  F: 52.8% M: 78.8%	Non-Parametric	C:41-49, D: 71-116, 163-209	hs-cTnI : WW except US Ver.3, Apr. 2021
LSI Medience PATHFAST cTnI-II	NP	NP	6.1%	NP	0.019	Overall n=333	Heparin-Na, heparin-Li or EDTA whole blood or plasma	Human cardiac I-T- C complex	Overall: 0.029	NP	Non-Parametric	C: 41-49, D: 71-116, 163-209	US: Ver.6, Dec.2022
QuidelOrtho/Triage/ Quidel TriageTrue hs- cTnI	0.4 ng/L (plasma) 0.5-0.8 ng/L (whole blood)	0.7 – 1.6 ng/L (plasma) 1.5-1.9 ng/L (whole blood)	5.0 – 5.9% at 21 ng/L (plasma) 5.9 – 6.5% at 22 ng/L (whole blood)	2.1 – 3.6 ng/L (plasma) 2.8 ng/L (whole blood)	4.4 – 8.4 ng/L (plasma) 5.8 – 6.2 ng/L (whole blood)	Overall n = 789  F: 391 M: 398	EDTA whole blood or plasma	NP	Overall: 20.5 ng/L  F: 14.4 ng/L M: 25.7 ng/L	Overall: ≥ 50%	NP	NP	hs-cTnI: WW except US April 2020
QuidelOrtho/Triage/ Quidel Triage Cardiac Panel	NP	0.050	NP	16.3% at 0.120	NP	Overall n = 323  F: 168 M: 155	EDTA whole blood or plasma	NP	NP	NP	NP	NP	USA, 2014, rev. D
QuidelOrtho/Triage/ Quidel Triage Profiler SOB Panel	NP	0.050	NP	16.3% at 0.120	NP	Overall n = 323  F: 168 M: 155	EDTA whole blood or plasma	NP	NP	NP	NP	NP	USA, 2014, rev. D
QuidelOrtho/Triage/ Quidel Triage Troponin I, Cardio2, Cardio3	0.002	0.01	NP	0.020	0.040	Overall n = 989	EDTA whole blood or plasma	NP	Overall: 0.020	Overall: 11.8%	Non-Parametric	NP	USA, 2014, rev. D
Radiometer AQT90 FLEX TnI	<0.005	<0.010	12.9%	NP	NP	Overall n = 231  F: 106 M: 125	EDTA and Li- heparinized whole blood and plasma	NP	Overall: 0.023	NP	Non-Parametric	C: 41-49, 190-196 D:137-149	Denmark 201808Z
Radiometer AQT90 FLEX TnT	0.005	0.010	15.2%	0.030	0.026	Overall n = 260  F: 132 M: 128	EDTA and Li- heparinized whole blood and plasma	Roche 4th gen	Overall: 0.017	NP	Non-Parametric	C:125-131 D:136-147	Denmark 201808J

Company/Platform/ Assay	LoB (µg/L unless noted)	LoD (µg/L unless noted)	% CV at 99 <sup>th</sup> Percentile	Conc at 20% CV (µg/L unless noted)	Conc at 10% CV (µg/L unless noted)	Reference Population N, Ages, Sex	Specimen Type	NP	99 <sup>th</sup> Percentile (µg/L unless noted)	Percent Normals Measured ≥ LoD	Statistic Used to Calc 99 <sup>th</sup> Percentile	Epitopes Recognized by Antibodies	Country of Package Insert: Version Date
Response Biomedical RAMP Troponin I	NP	0.03	20.0%	0.10	0.21	Overall n =180	Only EDTA whole blood	NP	Overall: <0.10	NP	NP	NP	Eu IFU 90012-1.2
Roche CARDIAC POC Troponin T – Roche cobas h 232	NP	0.04	NP	0.04 – 2.0	9.3% between 0.04 –0.2	Overall n = 302	Heparinized whole blood	NP	NP	NP	NP	D: 125-131 C: 136-147	EU, 2016-05
Siemens ATELLICA VTLi hs-cTnI <sup>#</sup>	0.55 ng/L	1.2 ng/L (plasma)  1.6 ng/L (whole blood)	6.5% (plasma)  6.1% (whole blood) at 22.9 ng/L	2.1 ng/L (plasma)  3.7 ng/L (whole blood)	6.7 ng/L (plasma)  8.9 ng/L (whole blood)	Overall n = 694  F: 331 M: 363	Li Hep whole blood and plasma, capillary blood	NP	Overall: 22.9 ng/L  F: 18.5 ng/L M: 27.1 ng/L	Overall: 84%  F: 80% M: 87%	Non-Parametric	C: 41-49 D: 2 Abs in range 20-100 cTnC Ab	CE, 2021-03
Siemens Stratus CS Acute Care cTnI test pack <sup>#</sup>	<0.03	NP	8.2%	0.03	0.06	Overall n=101	Whole blood (Li or NP heparin) or plasma Li or Na heparin	NP	Overall: 0.07	NP	NP	C: 27-32 D: 41-56	CE, 2008-04

LoB, limit of blank; LoD, limit of detection, NP, not provided; C, capture antibody; D, detection antibody; M, male, F, female; Conc, concentration; WW worldwide. All data have been listed as provided by the manufacturer, except assays in blue were abstracted from package insert due to lack of correspondence from manufacturer. <sup>#</sup> Disclaimer: The products/features mentioned herein are not commercially available in all countries. Their future availability cannot be guaranteed.